

AN INTRODUCTION TO OIM ANALYSIS™

Raising the
standard for
EBSD software

The most
powerful,
flexible, and
easy to use
tool for the
visualization
and analysis
of EBSD
mapping data



OIM Analysis™ - The Standard for EBSD Data Analysis

Why OIM Analysis™?


Orientation Imaging Microscopy (OIM) Analysis™ is a technique based on the automated collection and analysis of Electron Backscatter Diffraction (EBSD) patterns. This mapping data provides information on the orientation, phase distribution, grain size and shape, grain boundary structure, and local deformation of crystallographic microstructures. OIM Analysis™ provides virtually unlimited potential for interrogating the wealth of information contained in OIM scans.

OIM Maps

OIM Maps provide colorful and meaningful graphics, which give a useful insight into material microstructures.

The features of OIM Maps include:

- Greyscale (> 30 options), color (> 40 options) or combined maps, including orientation, phase, and EBSD image quality maps
- Overlay of misorientations, grain boundaries, phase boundaries, twin boundaries, CSL boundaries, and reconstructed boundaries
- Customization of map colors and parameters



OIM Charts

OIM Charts show numerical representations of the measured microstructure with statistical information for easy comparisons.

The features of OIM Charts include:

- More than 50 different chart options including grain size, misorientation distribution, and local misorientations analysis
- Correlated plots of measured parameters
- Customization of chart display

OIM Plots

OIM Plots offer visualization of measured orientation and misorientation distributions.

The features of OIM Plots include:

- Orientations are displayed in pole figures, inverse pole figures, Euler space, and Rodrigues space
- Misorientations are displayed in axis/angle space, Rodrigues space, and misorientation axis in pole figures
- Discrete, contoured or density plots


OIM Analysis™ - The Standard for EBSD Data Analysis

QuickGen Toolbar

The QuickGen Toolbar is an intuitive starting point for analysis offering quick access to commonly used maps, charts, and plots.

The features of the QuickGen Toolbar include:

- Quick orientation, image quality, grain, phase, and grain boundary maps
- Customized user-defined outputs
- Partition and dataset level one-button analysis



OIM Templates

OIM Templates offer an easy-to-use tool for performing repetitive, customized analysis of OIM data.

The features of OIM Templates include:


- Consistent reproduction of analysis and display parameters
- Templates for individual map/chart/plot, partition, and dataset levels
- Easy sharing of templates between users for improved collaboration

Data Partitioning


Data Partitioning enables the splitting of OIM data into different microstructural subsets so that the different portions can be measured separately.

The features of Data Partitioning include:

- Partitioning by point or by grain based properties
- Boolean logic to combine partition definitions
- Partition definitions can be saved as a template for repeated use



OIM Analysis™ - Powerful, Flexible and Easy to Use




OIM Project Tree


OIM Project Tree is an intuitive data management structure for novice and advanced users.

The features of OIM Project Tree include:

- Visual representation of current data analysis
- Cut and paste analysis within the project tree
- Creation and editing of maps, charts, plots and templates
- Access to data import, export (including HDF5 support), clean up, rotation, and other functions in the selected branch
- Both data and analysis can be saved together in the project
- OIM Data Miner for data management



Elastic stiffness map from additively manufactured titanium



Pole figure from textured aluminum


OIM Analysis™ - Powerful, Flexible and Easy to Use

Data Highlighting


Data Highlighting enables the analyst to identify specific data ranges and the following microstructural features on other maps and plots.

The features of Data Highlighting include:


- Highlighting from charts
- Highlighting from plots
- Highlighting from maps
- Adjustable color gradients for customized visualization
- Creation of partitions based on highlighted and non-highlighted regions for further analysis



— Microstructure
— Crack Path



Investigating cracking behavior in pipeline steels



Visualizing grain size distribution in aluminum

Interactive Analysis

Interactive Analysis allows the analyst to conduct a manual interrogation of data so that the information of interest can be extracted.

The features of Interactive Analysis include:

- Selection of specific points
- Selection of specific grains
- Selection of specific orientations
- Selection of specific regions of interest
- Measurement of point-to-point misorientations
- Measurement of triplet junction misorientations
- Measurement of vector point-to-point and point-to-origin data
- Demonstration of unit cell orientation and crystallographic plane geometry
- Ability to save and/or export interactive data


OIM Analysis™ - EBSD Results with Confidence

OIM Analysis™ EBSD Indexing

Even the best analysis tools may generate initial low-quality data. However, OIM Analysis™ has the resources to evaluate data quality and to improve it using EBSD pattern reindexing.

The features of EBSD Pattern Reindexing include:

- Reindexing data by individual points, by defined partitions, by phase or all data
- Reindexing using either saved band positions (Hough peaks) or using saved EBSD patterns
- Optimization of EBSD pattern band detection and indexing parameters for each phase
- Ability to find and add phases not included in the initial data collection
- Creation of EBSD signal background from individual phases and application of background processing optimized for specific materials
- Reindexing of multiple datasets using batch processing including Chl-Scan™ and NPAR™



NPAR™

NPAR™ is a patented approach for improving EBSD pattern indexing and mapping through local spatial pattern averaging.

The features of NPAR™ include:


- Improved EBSD pattern signal to noise ratio
- EBSD pattern acquisition at lower beam currents and/or faster camera speeds
- Higher EBSD mapping indexing success rates

Chl-Scan™

Chl-Scan™ offers a patented approach to improving phase differentiation and mapping using simultaneously collected EDS and EBSD data.

The features of Chl-Scan™ include:

- Automatic or manual chemical phase classification
- Selection of elements used for differentiation
- Use of PRIAS™ imaging signals as differentiation input




OIM Analysis™ - EBSD Results with Confidence

OIM Batch Processor

OIM Batch Processor enables the preparation and analysis of multiple datasets easily and consistently in a single action.

The features of OIM Batch Processor include:

- Batch analysis using user-defined OIM templates
- Alignment of correlated datasets using orientation information
- Ideal tool for analysis of in-situ data (heating stage, deformation stage, 3D)
- Batch cropping, data cleanup, and orientation rotations available
- Ability to save processed data as OIM Project for subsequent analysis





OIM 3D Visualization

OIM Analysis™ extends comprehensive EBSD characterization into 3D EBSD.

The features of OIM 3D Visualization include:

- Visualization of 3D EBSD data using OIM Map templates
- Orientations, EBSD image quality, phase, EDS elements, and local orientations
- Calculation of grains in 3D based on measured orientation information
- Calculation of grain size, grain average orientation, and grain orientation spread
- Transparency levels to visualize features within the microstructure
- Interactive selection and display of grains within the microstructure




OIM Applications

OIM has been used worldwide in thousands of publications and a wide range of applications.

Examples of OIM applications include:

- Improving the strength and ductility of lightweight automotive alloys
- Understanding the fatigue behavior of aerospace titanium alloys
- Tracing the evolution of trilobite eyes
- Increasing the efficiency of solar cells
- Increasing creep resistance in nuclear reactor materials



OIM Analysis™ Training

EDAX provides many options to assist users in learning how to use OIM Analysis™, including:

- Comprehensive manual
- Step-by-step tutorials
- How-to videos
- Training courses
- On-site training
- 'Ask the Expert' web sessions
- Topical webinars - live and on-demand

Summary

OIM Analysis™ is more than the sum of the functionality and specifications described in this brochure. In practice it is a powerful characterization and analysis tool, which enables researchers to better understand and evaluate the microstructure of materials.

To learn more about OIM Analysis™ or to see a demonstration, please visit our web site at www.edax.com or email us at info.edax@ametek.com.